

# HENSON 9000

## Key product features

<b>Fast, new, accurate algorithms</b>	New ZATA threshold algorithm conforms to Goldmann standard Faster threshold testing using prior data for greater accuracy Fast multiple and single stimuli suprathreshold program
<b>Data collaboration facilities</b>	Fully networkable through control device PDF format printouts available
<b>Standardisation</b>	Fully compliant DVLA estermann test Fully compatible with Henson 8000 Conforms with Goldmann standard
<b>Innovative technology</b>	Degradation-free light source Solid-state technology Low maintenance costs Numerous PC configuration options now available Optional motorised chinrest controlled by PC



## Technical specification

<b>Target type (Stimulator screen)</b>	Bowl	
<b>Target distance</b>	25 (cm)	
<b>Maximum stimulus illumination</b>	10,000 (Stimulus intensity) asb	
<b>Background illumination</b>	31.5 (asb)	
<b>Stimulus source</b>	LEDs white	
<b>Stimulus size</b>	Goldmann III	
<b>Presentation time</b>	Stimulus flash time (msec)	200
	Minimum inter-stimulus delay (sec)	0.5
	Patient response time	Adaptive or fixed
<b>Fixation</b>	Fixation target	Single and 4 point
	Heiji-Krakau	Yes
	Video eye monitor	Yes
<b>Test programs</b>	ZATA Standard threshold central	10-2, 30/24-2
	ZATA Fast threshold central	10-2, 30/24-2
	Suprathreshold single stimulus	1-3 Level. Can manually add test locations
	Suprathreshold multiple stimulus	1-3 Level. Can manually add test locations
	Estermann (DVLA driving test)	Binocular
<b>Unit dimensions</b>	W x D x H (mm)	
<b>Unit weight (kg)</b>	13.5kg	
<b>Input voltage</b>	85 - 263	
<b>Chinrest</b>	Yes	
<b>Headrest</b>	Yes	
<b>Database</b>	MS Windows™ compatible, networkable dependant on PC/laptop	
<b>Optional printer</b>	MS Windows™ compatible, dependant on PC/laptop	
<b>Control device</b>	External PC/laptop running Microsoft® Windows® 7 or 8.X Professional (32 or 64 bit)	



**Europe, Middle East and Africa**  
**Elektron Technology**  
Broers Building,  
JJ Thomson Avenue,  
Cambridge CB3 0FA UK  
T +44 (0)1223 371000  
W [elektron-healthcare.com](mailto:elektron-healthcare.com)  
E [info@elektron-healthcare.com](mailto:info@elektron-healthcare.com)

For more information about  
Elektron Technology go to:  
[www.elektron-technology.com](http://www.elektron-technology.com)

**The Americas**  
**Elektron Technology**  
5856 Corporate Ave, Suite 220,  
Cypress, CA 90630  
T +1-760-343-3650  
W [elektron-healthcare.com](mailto:elektron-healthcare.com)  
E [info@elektron-healthcare.com](mailto:info@elektron-healthcare.com)

**Asia Pacific**  
**Elektron Technology**  
3A International Business Park  
#08-02 ICON @ IBP, Tower A,  
Singapore 609935  
T +65 6933 7193  
W [elektron-healthcare.com](mailto:elektron-healthcare.com)  
E [info@elektron-healthcare.com](mailto:info@elektron-healthcare.com)



**Henson 9000.**  
Sensitivity. Specificity. Speed.  
Now you can have all three.

[www.elektron-healthcare.com](http://www.elektron-healthcare.com)



Henson. Screening and threshold testing with differences you need to see.

We are aware of the excellent technology available to ophthalmic and optometric professionals.

We are also aware of the frustrations and shortcomings of using this technology, in the context of practical and pragmatic patient care.

So we have created the **Henson 9000** that is designed to offer you a credible choice by offering you more of what you and your patients need. Without trade-offs.


**Small footprint, ease of use, performance and speed of operation have all influenced the new Henson 9000 design.**

The **Henson 9000** combines fast screening with the latest and fastest threshold test (ZATA) plus standard printouts that match those used in most Hospital Eye Services – see Fig. 1.

The modern and ergonomic design offers a headrest plus an optional motorised chin rest and a reduced footprint compared to previous models. There are multiple PC configuration options, enabling practices of any size to flexibly accommodate the unit.

When thinking about updating your perimeter, you really need to consider the major advantages the **Henson 9000** offers.

Screening: Henson introduces advantages impossible before.



(Fig 1) Test report printout

We know that when screening for disease, the two statistics determining performance are sensitivity and specificity.


We also know that a screening test needs to be conducted as quickly as possible, for the benefit of your practice and your patient.

Until now, these factors have been interdependent. For example, quicker testing has been achieved but at the expense of accuracy. Or increased accuracy has limited speed.

**High sensitivity?** We have scientific proof that you don't need large numbers of stimuli to have a test regarded as sensitive when testing for early visual field loss.


To access this paper please visit the *Henson Perimeter* section of [www.elektron-healthcare.com](http://www.elektron-healthcare.com)

**High specificity?** Our technology drives down false positive rates and allows clinicians to differentiate between random misses and glaucomatous defects. It enables you to repeatedly test a missed location.



**High speed?** The Henson screening test, in most cases, takes less than one minute per eye.

For a full in-depth explanation of our screening tests please visit the *Henson 9000* section of [www.elektron-healthcare.com](http://www.elektron-healthcare.com)



\*Not all tablets are suitable for use with the Henson 9000

Threshold testing: accuracy and speed are no longer mutually exclusive.

**ZATA. Different. And better.** For years, the SITA algorithm has been the engine underpinning threshold testing. Although SITA undoubtedly produces accurate test data, the way it works means testing is slow, which is undesirable for both clinician and patient. Attempts have been made to make SITA faster but there was a trade-off in accuracy.

ZATA – Zippy Adaptive Threshold Algorithm – was developed by Professor David Henson to overcome this apparently insurmountable issue. ZATA, although based on similar Bayesian methodology to SITA, manages to produce accurate results in a fraction of the time.

ZATA speeds up threshold testing by, where possible, using prior data. If there is some pre-existing visual field defect, this knowledge is taken advantage of and built upon. Rather than a test starting from scratch.

Secondly, ZATA doesn't just use single terminating criteria. It varies it to give more accurate thresholds at damaged and neighbouring locations.

This accelerates progress, particularly in patients with no visual field loss

Finally, ZATA uses looser determining criteria in severely damaged locations (<10dB). It doesn't even attempt to accurately measure these thresholds because it is impossible to do this with any algorithm.

For each and all of these reasons, the Henson 9000 with ZATA has significant advantages over other threshold testing devices.

For a full in-depth explanation of our ZATA tests, please visit the *Henson 9000* section of [www.elektron-healthcare.com](http://www.elektron-healthcare.com)



Facility	Benefit	Humphrey	Octopus
Uses Bayesian algorithm, ZATA	Faster, with no loss of accuracy	SITA	No Bayesian-based algorithm but alternative with loss of spatial resolution
Uses prior data	Faster, especially in patients with visual field loss. More efficient	Always starts from the same values (normal)	Not available
Varies the terminating criteria to give more accurate thresholds at damaged and neighbouring locations	Faster, especially in patients with no visual field loss	Uses the same terminating criteria for all locations	Not available
Uses looser terminating criteria in severely damaged locations (<10dB)	Faster, especially in patients who have severely damaged regions (common in hospital populations). Also reduces the number of non-seen presentations	Uses the same terminating criteria for all locations	Not available
30/24–2 stimulus patterns within single test. Operator can extend the 24–2 to a 30–2 test during or at the end of an examination	More efficient as patient does not have to repeat the 24-2 to do a 30-2 test	Has both 24–2 and 30–2 tests but not possible to extend from 24–2 to 30–2. If, at the end of a 24–2 test it is decided that a 30–2 is needed then you have to start again, repeating all the 24–2 test locations	Has both 24–2 and 30–2 tests but not possible to extend from 24–2 to 30–2. If, at the end of a 24–2 test it is decided that a 30–2 is needed then you have to start again, repeating all the 24–2 test locations
At end of test can switch view to threshold, grayscale, or defect values	All three views of the results can be seen on screen without having to print	Can only view threshold values on screen  Other indices are only available on the printout	

(Fig 3) Comparative chart for threshold tests