# Instrumententisch Instrument Table

ak 150 DL

Bedienungsanleitung  
User Manual

## INDEX

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
# User Manual

## INDEX

|           |   |           |
|-----------|---|-----------|
| <b>1</b>  | <b>GENERAL</b> .....  | <b>3</b>  |
| 1.1       | COPYRIGHT .....   | 3         |
| 1.2       | WARNING AND CAUTIONING SIGNS.....   | 3         |
| 1.3       | DISCLAIMER .....  | 3         |
| 1.4       | MANUFACTURER .....  | 3         |
| 1.5       | CE MARK, APPLICABLE RULES AND REGULATIONS .....                                   | 4         |
| 1.6       | USER INSTRUCTIONS FOR SAFE OPERATION AND OPERATOR QUALIFICATION REQUIREMENTS..... | 4         |
| <b>2</b>  | <b>USEFUL LIFE AND WARRANTY CONDITIONS</b> .....                                  | <b>5</b>  |
| <b>3</b>  | <b>SCOPE OF DELIVERY</b> .....  | <b>5</b>  |
| <b>4</b>  | <b>INTENDED USE</b> .....   | <b>5</b>  |
| <b>5</b>  | <b>MOUNTING AND COMMISSIONING</b> .....   | <b>6</b>  |
| <b>6</b>  | <b>ELECTRICAL CONNECTIONS</b> .....   | <b>6</b>  |
| <b>7</b>  | <b>PRODUCT DESCRIPTION AND CONTROLS</b> .....                                     | <b>7</b>  |
| 7.1       | POWER CONNECTION AND CONTROL UNIT.....  | 7         |
| 7.2       | CASTORS (LOCKABLE).....   | 7         |
| 7.3       | UP – DOWN CONTROL BUTTONS.....  | 8         |
| 7.4       | WARNING AND CAUTIONING SIGNS ON THE INSTRUMENT TABLE .....                        | 8         |
| <b>8</b>  | <b>OPERATION OF THE INSTRUMENT TABLE</b> .....                                    | <b>10</b> |
| 8.1       | DUTY CYCLE (CONTINUOUS OPERATION) OF ELECTRICAL MOTORS.....                       | 10        |
| 8.2       | ELECTRIC CONTROLS .....   | 10        |
| 8.2.1     | <i>Adjustment of table height</i> .....   | 10        |
| 8.3       | WARNING SIGNS WITH REGARD TO THE UP – DOWN MOVEMENT OF TABLE TOP .....            | 11        |
| <b>9</b>  | <b>ACCESSORIES</b> .....  | <b>12</b> |
| 9.1       | MOUNTING OF KEY BOARD DRAWER .....  | 12        |
| 9.2       | IEC SOCKET (2 PIECES) NOT RETROFITTABLE .....                                     | 12        |
| <b>10</b> | <b>CLEANING AND CONTAMINATION PROTECTION</b> .....                                | <b>13</b> |
| <b>11</b> | <b>MAINTENANCE</b> .....  | <b>13</b> |
| <b>12</b> | <b>DISPOSAL OF DEVICE</b> .....   | <b>13</b> |
| <b>13</b> | <b>TECHNICAL DATA</b> .....   | <b>14</b> |
| <b>14</b> | <b>ENVIRONMENTAL CONDITIONS</b> .....   | <b>14</b> |
| <b>15</b> | <b>TROUBLE SHOOTING</b> .....   | <b>15</b> |

## 1 General

To safely operate and control the medical device please carefully read this user manual before operating. Please pay attention to all cautionary statements pertaining to the safe operation of the medical device. Please keep the manual on file at all times while the product is in service.

The manufacturer reserves the right to make technical changes or modifications in accordance with advanced technical developments. This manual is not subject to an updating service.

### 1.1 Copyright

© Distributing, copying or any commercial use of this manual is strictly prohibited unless expressly permitted in writing by the manufacturer. The manufacturer is entitled to compensation for any violation of this right.

The manufacturer reserves any and all rights in case a patent is granted or the industrial design is registered.

### 1.2 Warning and cautioning signs

|   |         |  |
|---|---------|--|
|   | Warning | Indicates a hazardous situation which may result in fatal or serious bodily injury if the appropriate safety precautions are not heeded.   |
|  | Caution | <p>Indicates a hazardous situation which may result in minor injury if the appropriate safety precautions are not heeded.</p> <p>Caution Property damage<br/>Indicates possible device damage if the appropriate safety precautions are not heeded</p> |
|  |         | This pictograph on the label indicates Device typ B according to EN 60601-1  |

### 1.3 Disclaimer

|   |         |   |
|---|---------|---|
|  | Warning | This device may only be set up, operated and used for the specified purpose and according to national regulations, consistent with the applicable industry standards and occupational safety and accident prevention regulations. |
|---|---------|---|

This declaration shall be rendered invalid if changes are made to the product without the manufacturer's authorization.

### 1.4 Manufacturer

#### AKRUS GmbH & Co KG

Otto-Hahn-Straße 3

D-25337 ELMISHORN

int. +49 (0) 4121 791930

FAX int. +49 (0) 4121 791939

Email: [info@akrus.de](mailto:info@akrus.de)

Website: [www.akrus.de](http://www.akrus.de)

1.5 CE mark, applicable rules and regulations

The device complies with the following Regulations for machine products.

- 2006/42/EG
- RL 2011/65/ EU

The device meets the requirements of 2006/42/EG and the corresponding national requirements of the machinery directive.

User must observe the national legal requirements governing the prevention of accidents.

1.6 User instructions for safe operation and operator qualification requirements

**Before using this medical device, please carefully read and observe the safety instructions and recommendations listed in this manual**

|   |   |  |
|---|---|--|
|  | <p>Caution<br/>Risk of operating errors</p> | <p>The device may only be installed, operated, used and maintained by persons who have been properly trained or who have the required knowledge and experience to do so. Please also adhere to the national qualification guidelines applicable in your country.</p> |
|---|---|--|

- Proper handling and use of the medical device is essential to its safe operation.
- Operation of the medical device only in accordance with the intended use.
- Do not operate the medical device in explosive or hazardous environments. The medical device must not be operated, with the exception of very minute quantities, in the presence of combustible anesthetics, or volatile solvents like alcohol, benzene or the like.
- Do not set up the medical device in humid or damp rooms. Never expose the product to sprinkling, dripping or splashing water.
- Pay special attention to and carefully read any emphasized safety advice or information in this user manual.

## 2 Useful life and warranty conditions

The useful product life is expected to be 8 years. Both, the useful life time and the warranty offered by manufacturer are conditional upon the adherence to the conditions outlined below.

This declaration shall be rendered invalid if changes are made to the product without the manufacturer's authorization.

|   |   |   |
|---|---|---|
|    | <p>Caution<br/>Risk of operating errors</p> | <p>The device may only be installed, operated, used and maintained by persons who have been properly trained or who have the required knowledge and experience to do so. Please also adhere to the national qualification guidelines applicable in your country.</p>  |
|    | <p>Warning</p>                              | <p>The development, production and maintenance of these instrument tables, together with associated risks, are based on an expected service life of eight years, provided the device is serviced at the specified intervals.<br/>Modifications to the product or failure to follow the manufacturer's instructions may substantially reduce the expected service life and significantly increase the risks associated with the use of this device. It is the responsibility of the institution operating this product to follow the manufacturer's instructions and to decide on the risk/benefit ratio upon expiration of the expected service life or maintenance and inspection intervals specified by the manufacturer.</p> |
|  | <p>Warning<br/>General Hazards</p>          | <p>This device may only be set up, operated and used for the specified purpose and according to national regulations, consistent with the applicable industry standards and occupational safety and accident prevention regulations. Further notes on classification are to be found in chapter <i>Technical data</i>.</p>  |

## 3 Scope of delivery

The instrument table scope of delivery comprises

| Position                   | Qty |
|----------------------------|-----|
| Instrument table ak 150 DL | 1   |
| Power cord                 | 1   |
| Product documentation      | 1   |

## 4 Intended use

- The instrument table ak 150 DL is designed to carry ophthalmic medical devices and systems up to a max. work load of 100 kg.
- The product is suitable to be employed in the immediate proximity of patients.

## 5 Mounting and commissioning

The medical device will be delivered either completely assembled or folded down as an assembly kit. If not expressly mentioned in this user manual, no additional assembly or calibration work is required to operate the fully assembled medical device.

In case of the assembly kit, please follow the assembly instruction, which are part of the shipment

## 6 Electrical connections


The location of the power input of the device is described in this user manual. To connect the medical device to the power supply use only the power cable included in this delivery.

|   |   |  |
|---|---|--|
|    | <p>Warning<br/>General Hazard</p>         | <p>Do only store or operate the medical device in compliance with the environmental conditions as outlined in section 14.<br/>The device should be set up so that the power cable can be disconnected from the power supply quickly and easily without any tools.</p>  |
|    | <p>Warning<br/>Risk of electric shock</p> | <p>Additional extension cables or portable multiple sockets may not be used.<br/>The electrical installation must conform to IEC 60364-7-710 or the applicable national regulations. This includes the existence of a ground fault circuit interrupter (GFCI).<br/>To avert the risk of an electric shock this device may only be connected to a power supply with a protective earth conductor.<br/>Ensure that the power supply plug is suitable and certified for the local connection. If the supplied power cable must be replaced, at least the following specifications must be adhered to:</p> <ul style="list-style-type: none"> <li>• Protective ground conductor resistance of the power cable maximum 0.1 Ohm</li> <li>• Local certification of the power cable for connection to medical devices</li> <li>• Device plug C19 conforming to IEC 60320</li> <li>• Cross-section at least 1.0 mm<sup>2</sup></li> </ul> |
|  | <p>Warning<br/>Fire Hazard</p>            | <p>The device is not suitable for operation in explosion risk areas (e.g. combustible mixture of anesthetic, cleaning or disinfecting agents with air, oxygen or nitrous oxide).<br/>The electrical installation must conform to IEC 60364-7-710.<br/>The data concerning power consumption on the type plate must be taken into account in the selection of over current protection.</p>  |

## 7 Product description and controls

### 7.1 Power connection and control unit

The instrument table has a wide range power supply unit inside the table top adapter

### 7.2 Castors (lockable)

The instrument table can be moved on castors. To safely lock the table in position, make sure the castors are locked.

### 7.3 Up – Down control buttons



The rocker switch to move the table up or down are inside of the table top holder. The rocker switch is marked with up – down symbols.

1. UP / DOWN switch

### 7.4 Warning and cautioning signs on the instrument table

|          |  |   |
|----------|--|---|
| <p>1</p> |   | <p> Manufacturer</p> <p> Date of manufacture</p> <p> EU- Conformity Symbol</p> <p> AC current</p> <p><b>IP21</b> Ingress protection rating of casing (protected against foreign solid bodies <math>\geq</math> 12.5 mm in diameter and water drips from above)</p> <p> Use permitted only indoors</p> <p> Do not dispose of in domestic waste</p> <p><b>REF</b> Catalogue number / part number</p> <p><b>SN</b> Serial No.</p> <p> Danger to patient, operator</p> |
| <p>2</p> |   | <p>Information sign "Observe user manual"</p>   |

## 8 Operation of the instrument table

### 8.1 Duty cycle (continuous operation) of electrical motors

The electrical motors are designed for a duty cycle of approximately 2 minute at full load. Continuous operations exceeding that time period may cause overheat conditions and permanent damage to the transformer. (refer to section 8.3)

### 8.2 Electric controls

The table height can be adjusted (UP-DOWN) by pushing the corresponding control buttons. (refer to section 7.3)

#### 8.2.1 Adjustment of table height

The table top is height adjustable by two electrical lifts. Set table top to line up the instruments to match the body size of the patient. Upon release of the button, the table top will stop traveling at that position. Do not push both control buttons at the same time as this may cause damage to the electric lift.

### 8.3 Warning signs with regard to the Up – DOWN movement of table top


In case of very tall wheelchair patients (> 1,90m) make sure that, during downward movements of the table, sufficient space between patients thighs and the moving table parts is available. Danger of pinching. After an operating period (continuous duty cycle) of 2 minute a cooling off period of 20 minutes is required.

## 9 Accessories

### 9.1 Mounting of key board drawer

To mount key board drawer to the table, follow the assembly instructions below.  
Fig. key board drawer

### 9.2 IEC socket (2 pieces) not retrofittable

## 10 Cleaning and contamination protection

- Contamination protection  
In order to protect the product against contamination by liquids, body liquids or other unwanted substances it is recommended to, while in use, cover the product with a non sterile disposable impermeable cover sheet.
- Surfaces are resistant to all commonly used disinfectants as listed in the D G H M (Deutsche Gesellschaft fuer Hygiene und Mikrobiologie) list. [www.dghm.org](http://www.dghm.org)
- Do not attempt to sterilize the medical device.
- Clean only the external surfaces of the equipment using a damp cloth. In case of heavy stains do **not** use abrasive or aggressive materials other than commonly used cleansers or detergents. Do not use ether, acetone or concentrated acids.
- No water or cleaning detergents may enter the inside of the equipment.
- To clean soiled upholstery surfaces, use a damp (not wet) cloth and/or commonly used upholstery spray.
- Before resuming operation, allow sufficient time to dry.

## 11 Maintenance

The medical device is maintenance-free to users and operators.

The product is subject to the regular service intervals. Refer to the current issue of manufacturers "maintenance record".

## 12 Disposal of device

In accordance with applicable EU guidelines and national regulations at the time at which the product was brought onto the market, the product specified on the consignment note is not to be disposed of via the domestic waste disposal system or communal waste disposal facilities.

For further information on disposing of this product please contact your local dealer or the manufacturer or its legal successor company. Please read the latest information provided on the manufacturer's website.

If the product or its components are resold, the seller must inform the buyer that the product must be disposed of in accordance with the currently applicable national regulations.

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
### 14 Environmental Conditions

| <b>Ambient conditions for intended use</b>                                |                             |
|---|-----------------------------|
| Ambient air temperature   | +10° C - +40° C             |
| Rel. Humidity max   | 50 % non-condensing         |
| Ambient air pressure  | 700 – 1060 hPa              |
| <b>Ambient conditions for storage</b>                                     |                             |
| Ambient air temperature   | -10 °C to +55 °C            |
| Rel. Humidity   | 10 % to 95 % non-condensing |
| <b>Ambient conditions for storage and transport in original packaging</b> |                             |
| Ambient air temperature   | -40 °C to +70 °C            |
| Rel. Humidity   | 10 % to 95 % non-condensing |

## 15 Trouble shooting

| Fault/Malfunction     | Possible cause  | Corrective action   |
|-----------------------|---|---|
| No function available | Power cord disconnected/defective   | Check power cord  |
|                       | Up – Down buttons defective   | Replace Up – Down buttons, call service   |
|                       | Electric lift(s) or control unit defective  | Replace lift, by authorized service staff only  |
|                       | The lifting columns (control) must be re-initialized<br> | For initialisation, the button must be pressed for at least 10 seconds kept pressed during and the entire initialisation run. During the initialisation run, the Multilifts are traversed to the lowest position at half speed and not synchronously. A buzzer sounds five times at the lowest position - thereafter the drives are extracted to the highest position synchronously at normal speed. In the highest position, the buzzer sounds again and goes out after releasing the button |
| Mechanical damages    | External forces   | Call service  |